## **Ellis Mini-Fixateur extern**

## **Ellis External Mini Fixator**

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**Osteosynthesis: Fixator and Instruments**

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
Fixierdrähte/Fixator pins

**The new LINK Ellis External Mini Fixator has been developed for surgery of the:**

- Hand
- Wrist/distal radius
- Ulnar shaft
- Foot

**Indications:**


1. Fractures
2. Osteotomies
3. Arthrodeses
4. Non-union repair
5. Bone tumor surgery

**The major features:**


- Small size and low weight
- Direct pin-to-rod fixation
- Simple and strong construction
- Versatility
- Low cost

**The basic three components:**

1. Fixator pins
2. Threaded rods with one or more slots
3. Fixation nuts

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
**Bone lengthening:**

With a lengthening component for lengthening procedures at small bones.

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### Fixator construction:

The pins match precisely the slot of the rods. To set up the external fixator, the pins must have passed the slot in perpendicular direction. Fixation is achieved by tightening the nuts on the rod on both sides of the protruding pins.

### The major differences to other devices are:

1. **Direct pin-to-bone-fixation** without connecting device.
2. **The possibility to bend the pins.** This offers the advantage of oblique pin fixation of a fracture or osteotomy and provides initial stabilization.  
A comprehensive variety of external fixation configurations can be created for unilateral or bilateral rod application (transfixation).  
See illustrations on page 4.
3. **Lengthening components optional.**

### Two types of pins are used:

1. **Smooth pins.** Used whenever the pins have to be bent or in transfixation.
2. **Pins with a threaded tip.** Used in unilateral fixator application for the perpendicular positions to achieve sufficient purchase (grip) in the bone and to avoid pull-out.

### The external fixator is available in two system sizes:

- **1.5 mm** pins and rods for handsurgery.
- **1.9 mm** pins and rods for wrist, ulnar shaft and foot-surgery.

### Options:


Connectors, either straight or angled, to connect multiple rods.

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**General technique for bending the pins:**

1. If any oblique pin set-up is considered, the pin should be placed in the bone first.
2. After placement, hold the pin tightly with the flat nosed pliers to prevent tearing out.
3. Slide the bending tube over the pin, until it touches the pliers nose.
4. Bend the pin until it is perpendicular to the bone.
5. Slide a rod without nuts onto the pin as a drill guide for placing the additional pins.
6. Following wire introduction, the threaded rod is removed, completed with the necessary number of fixation nuts (2 ea./wire) and again applied. Finally the wires are fixed as described.

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**Fixator configuration models:**


- Fig.1A Phalangeal, metacarpal and metatarsal fractures and osteotomies. Ulnar shaft fractures.
- Fig.1B Transfixation in comminuted phalangeal fractures, tumour resections or in case of poor bone quality.
- Fig.1C Oblique fractures of phalanges, metacarpals and metatarsals. Trapeziometacarpal arthrodeses.
- Fig.1D Arthrodeses of first metatarso-phalangeal and metacarpo-phalangeal joint. Distal radial fractures and osteotomies.
- Fig.1E Transfixation with crossed pins for arthrodesis of interphalangeal joint.
- Fig.1F Example of fixation of a metaphyseal fracture.

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Trapeziometacarpal arthrodeses with external fixation offers the advantage of full wrist and thumb movement during treatment (without cast) which cannot be achieved with other types of osteosynthesis (Fig. 1, 2, 3).

### The fixator may be applied without additional immobilization in:

- Phalangeal fractures and osteotomies of the hand and foot
- Metacarpal fractures and osteotomies
- Phalangeal non-union repair
- Interphalangeal arthrodeses of fingers
- Trapeziometacarpal arthrodeses
- Distal radial osteotomies
- Open ulnar shaft fractures
- First metatarsal osteotomies (Hohmann)
- First metatarsophalangeal arthrodesis
- etc.

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**Fig. 1** Model of trapeziometacarpal arthrodesis, according to configuration model no.1C (page 4).

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**Fig. 2** External fixation in trapeziometacarpal arthrodesis after malunited first metacarpal osteotomy. Full wrist and thumb motion is ensured.

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**Fig. 3** After removal of the fixator: consolidation in a functional position of the thumb.

**Schuind, F.**

Small External Fixation Devices for the Hand and Wrist.  
*Clin Orthop* 293, 1993.

**Bamberger, H.**

Trapeziometacarpal Joint Arthrodesis: A Functional Evaluation.  
*J Hand Surg* 17A, 1992.

**Caputo, R.**

Power Staple Fixation in Trapeziometacarpal Arthrodesis.  
*J Hand Surg* 18A, 1993

**Ellis, M.L.**

Hands on the hand. The development of a new external mini fixator.  
*Read at the Annual Meeting of the Dutch Orthopedic Society, Utrecht, The Netherlands, October 8, 1994.*

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
- a) Painful fifth metatarsal head with clavus on the little toe due to the varus position.
- b) After wedge osteotomy with medial shift and external fixation with a three-pin half-external fixator.
- c) After removal of the fixator and consolidation of the osteotomy.

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**Lateral view.**  
Nice alignment of the toes.

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**Plantar view of the operated foot.** Full weight bearing was possible.

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**Frontal view.** Smaller circumference of the foot at the level of the metatarsal heads.

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

Unstable fracture of the proximal phalanx of the fifth digit.

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
Anatomical position after healing. Six weeks after fixation.

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**Bone lengthening device**

The **LINK® Ellis External Mini Fixator** system is available with a bone lengthening component. This device has been developed to enable lengthening procedures of small tubular bones. It has been designed to use together with all other parts of the Ellis Fixator. A lengthening device for either 1.5 mm and 1.9 mm fixation pins is available.

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Indications for use include congenital or posttraumatic shortening of digits, metacarpals and metatarsals.

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
The bone lengthening technique is illustrated on a second metacarpal using 1.5 mm fixation pins.

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**Pin placement and osteotomy**

A small longitudinal incision is made to expose the metacarpal shaft. Four fixation pins with threaded tip are drilled parallel to each other, two proximal and two distal to the designated osteotomy site. In each pair, the pins should be placed as far from each other as possible. Then the osteotomy should be performed according to the surgeons preference.

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

Insertion of four parallel fixation pins with threaded tips. Then osteotomy of the second metacarpal.

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
### Mounting of the Ellis Fixator and the bone lengthening device

On the proximal and distal fixation pins two threaded rods of 32 mm are fixed with fixation nuts. Additional nuts to hold the lengthening component must be placed on the threaded rods in advance with the fixation nuts. As a next step the lengthening component has to be prepared. Lengthening components are available in lengths of 35, 45 and 55 mm. The necessary component length depends mainly on the desired bone length and the placement of fixation pins. First a numbered nut, used to define the individual steps of distraction, is screwed on the lengthening component up to the fixed coupling rod. Successively the adjustable coupling rod


is placed on the extension rod. Prongs of the fixed as well as the adjustable coupling rod are now placed over the fixator's threaded rods. The fixed coupling rod should be situated proximally. A standard fixation nut is then screwed on at the distal end of the lengthening component. Finally the adjustable coupling rod between numbered nut and standard fixation nut is fixed with measured compression force. The set-up is now stable.

### Bone lengthening

Bone lengthening is performed by loosening the distal nut of the hand turning the broad nut to the right (from the patient's point of view). The figures indicate the amount of lengthening. One full 360° turn of the nut produces 0.8 mm of

lengthening. One increment of the figure on the nut corresponds therefore with 0.13 mm. Lengthening of 1 cm requires 12.5 full rotations of the nut.

Lengthening can be achieved to any desired distance.

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**Fixator set and instruments, complete**  
on two storage trays, in container

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**Threaded rods for wrist, ulna and  
foot surgery**

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




**Lengthening component** mounted on two threaded rods with three fixation pins.

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03/04/2009

**Please note the following regarding the use of our implants:**

- 1. Choosing the right implant is extremely important.**  
The size and shape of the human bone determine the size and shape of the implant and also limit the load capacity. Implants are not designed to withstand unlimited physical stress. Demands should not exceed normal functional loads.
- 2. Correct handling of the implant is exceedingly important.**  
Under no circumstances should the shape of a finished implant be altered, as this shortens its life span. Our implants must not be combined with implants from other manufacturers.  
  
The instruments indicated in the Surgical Technique must be used to ensure safe implantation of the components.
- 3. Implants must not be reused.**  
Implants are supplied sterile and are intended for single use only. Used implants must not be reused.
- 4. After-treatment is also very important.**  
The patient must be informed of the limitations of the implant. The load capacity of an implant cannot compare with that of healthy bone!
- 5. Unless otherwise indicated, implants are supplied in sterile packaging.**  
Note the following conditions for storage of packaged implants:
  - Avoid extreme or sudden changes in temperature.
  - Sterile implants in their original, intact protective packaging may be stored in permanent buildings up until the „Use by“ date indicated on the packaging.
  - They must not be exposed to frost, damp and direct sunlight or mechanical damage.
  - Implants may be stored in their original packaging for up to 5 years after the date of manufacture. The “Use by” date is indicated on the product label.
  - Do not use an implant if the packaging is damaged.
- 6. Traceability is important.**  
Please use the documentation stickers provided to ensure traceability.
- 7. Further information** on the material composition is available on request from the manufacturer.


### **Follow the instructions for use!**

WALDEMAR LINK GmbH & Co. KG, Hamburg.

All content in this catalogue, including text, pictures and data, is copyright-protected. Every instance of use not permitted by the German Copyright Act is subject to our prior consent. In particular, this applies to the reproduction, editing, translation, saving, processing or passing on of content stored in databases or other electronic media and systems. The information in this catalogue is solely intended to describe the products and does not constitute a guarantee.

The Surgical Technique described has been written to the best of our knowledge and belief but it does not relieve the surgeon of his responsibility to duly consider the particularities of each individual case.

Unless otherwise indicated, all instruments are made of surgical stainless steel.

WALDEMAR LINK GmbH & Co. KG  
Barkhausenweg 10 · D-22339 Hamburg  
Postfach 63 05 52 · D-22315 Hamburg  
Tel +49 (0)40 5 39 95-0 · Fax +49 (0)40 5 38 69 29  
e-mail [info@linkhh.de](mailto:info@linkhh.de) · Internet [www.linkhh.de](http://www.linkhh.de)